

# Zoll M series SAED Battery Management Protocol

Follow this procedure each fortnight *and* following extended use of the Zoll ....

- 1. Firstly, remove the Battery (A) which is currently in the monitor.
- 2. Take the spare battery (B) from the Zoll pack and insert into monitor




**Rotation of batteries**  
Where Battery 'A' is charged one fortnight, Battery 'B' will be charged in the same way the next fortnight...



Load battery (A) into the *Power Base Charger* & Press the 'Test' button

**Hint...**  
Mark each Battery as "A" or "B" to distinguish between them



Leave the battery to test and charge. This takes ~ 8 hours

Place fully charged Battery in the rear pocket of the Zoll pack

Only when the 'Ready' light is illuminated Remove battery 'A'

Write the date on the back of the battery when this charge & test procedure is completed